assero H20P - The very best for moisture in oil.

assero H20P Portable Moisture in Oil Meter

For the very best results from your moisture in oil analysis, the assero H2OP has been developed in response to many comments from our customers in laboratories and in the field, to develop a portable, robust, easy to use and highly accurate moisture meter for the measurement of down to 1ppm H2O in insulating oils.

The microprocessor controlled assero H20P has a built-in battery and optional carry case , provide the versatility required by the laboratory and also the ease of use and portability required by the field personnel and carries a 5 year warranty.

Easy to Use

Simple to programme so that only a single button needs to be pressed for the analysis , everything else is automatic.

The unique hermetically sealed glassware design is designed to eliminate any outside atmospheric moisture getting into the cell, allowing the highest sensitivity and performance on the market. The "quick-connect" system hermetically seals the glassware without the use of expensive grease or Teflon sleeves allowing trouble free assembly, maintenance and disassembly.

Results & Records

The built-in high speed printer provides a hard copy of results and allows duplicate print outs if required. Results can also be downloaded via the Results Manager software package onto a pc spreadsheet format.

Spares & Accessories

The assero H2OP is supplied ready for operation. A comprehensive glassware pack including electrodes, Glass titration vessel, leads, syringe, printer paper, etc, is supplied as standard. A calibration certificate and 5 year warranty is also provided.

Reagents

Specially formulated to provide optimum performance with all makes and models of coulometric Karl Fischer titrators, our unique reagents provide improved sample miscibility and solubility. assero H2OP reage nts are packaged in single-shot screw top bottles containing 100ml anode reagent and 5ml easy-snap ampoules of cathode reagent. A standard carton contains 8 anode charges and 8 cathode charges. The reagent pack comprises 8 x 100ml anode reagent and 8 x 5ml cathode reagent.

assero H2OP Coulometric KF Moisture in Oil Meter \cdot Simple operation \cdot 10 user programmable methods \cdot 1ppm / 100% \cdot Results in ppm, mg/Kg, % & µg water $\neg \cdot$ Multi language display & print out $\neg \cdot$ Small footprint $\neg \cdot$ Integral high speed printer $\neg \cdot$ Fully portable $\neg \cdot$ Low drift cell design $\neg \cdot$ Automatically compensated errors (patented technique) \cdot Optional carry case $\neg \cdot$ Conforms to ASTM D1533, D 4928, D6304, IP386, IP438, API MPMS Chapter 10.9, IEC60814, ISO 10101-3, 10337 & 12937.







Technical Specification - assero H20P

Titration method	Coulometric Karl Fischer titration
Electrolysis control	Automatic Compensated Errors control system
End point detection	AC polarisation
End point indication	Visual display/print out/acoustic beep
Display	40 character alphanumeric backlit LCD
Measuring range (possible)	1µg – 100mg water
Measuring range (typical)	1µg – 10mg water
Moisture range	1 ppm – 100% water
Max. sensitivity	0.1 µg
Max. titration speed	2.0 mg per minute
Max. current	400 ma
Drift compensation	Automatically controlled
Start delay time	0 – 30 minutes, user selectable
End delay time	0 – 30 minutes, user selectable
Power supply	90-264V AC, 47-63 Hz. 12V DC car adapter/internal battery

Precision	10-100µg ±3µg, 100µg-1mg ±5µg, above 1mg ±0.5%
Calculation modes	Weight/weight, user programmable
	Weight/dilution ratio, user programmable
	Volume/density, user programmable
	Volume/volume, user programmable
Display format	μg, mg/kg, ppm, %
Print format	μg, mg/kg, ppm, %
Statistics	max, mean, min values upto 99 runs
Method storage	10 user programmable methods
USB output	USB Flash Drive
Sample ID number	User programmable
Printer	42 character high speed thermal printer
Stirrer speed	Microprocessor controlled
Dimensions & Weight	250 x 245 x 120 mm, 3kg
Calendar/clock	Analysis time & date print out
Battery low indicator	Display & print out indication